

CTDS 293 SFM Basic

CTDS 293 SFM Basic is a serum-free medium developed for HEK293 cells.

Specialty culture media help tailor the extracellular environment to a defined cell type or research workflow.

Appropriate formulation selection can simplify the development of consistent culture conditions.

CTDS SKU	CTDSK300
Pack size	1000 ml
Product category	Cell Culture Media
Product subcategory	Protein Expression Media
Detailed category	Transient Expression Media

PRODUCT BACKGROUND

CTDS 293 SFM is a serum-free medium developed for HEK293 cells. It is suitable for culture, scale-up or production-oriented workflows depending on the cell line, supporting high-density growth, stable performance and reduced serum-associated variability. It can be used for recombinant protein expression, vaccine process development, cell line screening or process optimization after validation. For research, process development and manufacturing support use; final supplementation should be optimized as per cell line and protocol.

APPLICATIONS

Used in cell culture applications

CHARACTERIZATION AND QUALITY INFORMATION

Assay / validation method	Flow cytometry marker characterization
---------------------------	--

HANDLING AND STORAGE

Shelf life	9 months
Shipping	Shipped on gel pack
Storage	Store at 2-8°C
Intended use	For Research, Laboratory and Industrial Use Only. Not for Diagnostic Use

PRODUCT DISCLAIMER

The data and recommendations are based on validation done in our laboratory. It is recommended for users to establish their own parameters of culture independent of our claims. These reagents are for research use only and not for in-vitro diagnostics or human use.

ADDITIONAL CLASSIFICATION

Product grouping	Serum-Free Media (SFM)
Product subgroup	Serum-Free Medium (Individual Components) For Transient Protein Expression

Not a lot-specific certificate of analysis.