

CTDS Gln Max

CTDS Gln Max is a stabilized L-alanyl-L-glutamine supplement designed to provide a durable glutamine source for cell culture media.

Culture supplements help adjust the nutrient environment without replacing the complete basal formulation.

Their use allows researchers to tailor culture conditions to the requirements of a particular cell system.

CTDS SKU	CTDSK900
Pack size	100 ml
Product category	Sera and Supplements
Product subcategory	Culture Supplements

PRODUCT BACKGROUND

CTDS Gln Max is a stabilized L-alanyl-L-glutamine supplement designed to provide a durable glutamine source for cell culture media. It helps reduce instability associated with conventional glutamine and supports improved culture consistency in serum-free, low-serum and standard mammalian culture systems. For research, process development and manufacturing support use; final supplementation should be optimized as per cell line and protocol.

APPLICATIONS

Used in cell culture applications

CHARACTERIZATION AND QUALITY INFORMATION

Assay / validation method	Flow cytometry marker characterization
---------------------------	--

HANDLING AND STORAGE

Shelf life	12 months
Shipping	Shipped on gel pack
Storage	Store at 2-8°C
Intended use	For Research, Laboratory and Industrial Use Only. Not for Diagnostic Use

PRODUCT DISCLAIMER

The data and recommendations are based on validation done in our laboratory. It is recommended for users to establish their own parameters of culture independent of our claims. These reagents are for research use only and not for in-vitro diagnostics or human use.

ADDITIONAL CLASSIFICATION

Product grouping	Serum-Free Media (SFM)
Product subgroup	Other Medias / Reagents

Not a lot-specific certificate of analysis.