

CTDS ITS

CTDS ITS is an insulin-transferrin-selenium supplement for serum-reduced and serum-free cell culture applications. Specialty culture media help tailor the extracellular environment to a defined cell type or research workflow. Appropriate formulation selection can simplify the development of consistent culture conditions.

CTDS SKU	CTDSK930
Pack size	10 ml
Product category	Sera and Supplements
Product subcategory	Culture Supplements

PRODUCT BACKGROUND

CTDS ITS is an insulin-transferrin-selenium supplement for serum-reduced and serum-free cell culture applications. It supports growth and maintenance of many mammalian cells by providing defined nutritional and trace element support while helping reduce reliance on serum. For research, process development and manufacturing support use; final supplementation should be optimized as per cell line and protocol.

APPLICATIONS

Used in cell culture applications

CHARACTERIZATION AND QUALITY INFORMATION

Assay / validation method	Flow cytometry marker characterization
---------------------------	--

HANDLING AND STORAGE

Shelf life	12 months
Shipping	Shipped on gel pack
Storage	Store at 2-8°C
Intended use	For Research, Laboratory and Industrial Use Only. Not for Diagnostic Use

PRODUCT DISCLAIMER

The data and recommendations are based on validation done in our laboratory. It is recommended for users to establish their own parameters of culture independent of our claims. These reagents are for research use only and not for in-vitro diagnostics or human use.

ADDITIONAL CLASSIFICATION

Product grouping	Serum-Free Media (SFM)
Product subgroup	Other Medias / Reagents

Not a lot-specific certificate of analysis.